

## Metal Chelating Wipes (IPA/Water/Chelating Agent Lint-Free) 100 wipes, 9"x6.1" each

### Product Highlights

Metal Chelating Wipes  
Helps remove metal from surfaces and skin by forming water-soluble metal-citrate compounds with metals present to allow them to be wiped off

**RoHS 3 and REACH compliant**

### Specifications

#### Formula:

70% IPA / 30% Deionized Water + Chelating Agent

#### Packaging:

Cannister (100 Lint-Free Saturated Wipes, 9" x 6.1" (23cm x 15.5cm)

#### Net Liquid Weight Saturated into Wipes:

550g

### Shelf Life

Refrigerated >24 months, Unrefrigerated >24 months

### Storage and Handling

Store refrigerated or at room temperature 3-25°C (37-77°F). Do not freeze.

### Transportation

Flammable. Ground Shipping Only. Shipping below 0°C (32°F) or above 25°C (77°F) for normal transit times by ground or air will not impact this product's stated shelf life.

### Conforms to the following Industry Standards:

RoHS 3 Directive (EU) 2015/863:

Yes

